




# ssR Amino Linkers



## ssR Amino Linkers: ssH & ssMe Amidites



ssH Amino linker; CLP-1132




ssMe Amino linker; CLP-1131

## Advantages and Key Features:

ssR amino linkers with an aminoethyl carbamate main linkage and a side chain residue have several following advantages over conventional amino linkers.


- Strong hydrophobicity of MMTr group assists in easy purification of the amino modified oligonucleotides.
- High solubility in acetonitrile and standard oligonucleotides synthesis protocols can be applied and it allows measurement of the coupling efficiency through colorimetric MMT-assays.
- The aminoethyl carbamate structure facilitates deprotection of MMTr group under very mild acidic conditions compared to the standard aliphatic amine and this feature avoids the depurination side reaction.
- ssR amino modifications have displayed very high labeling efficiency with an active esters, isothiocyanate compared. (1) with Biotin-NHS/phosphate buffer (pH8), 77-79% yield for ssH and 82-84% yield for ssMe, with FITC/phosphate buffer (pH8), 57-59% yield for ssH and 63-65% yield for ssMe reporter group at the 5' end of oligonucleotides.
- **Pegylation** with polyethylene glycol: 40K PEG-NHS\* with ssH-linker yield is 70-75%.
- **Pegylation** 60% for ssH and 65.6% for ssMe amino linkers.
- **Cholesterol** attachment coupling efficiency is 70-80%.
- High-throughput purification due to increased conjugation efficiency.

## Comparative illustration of pegylation with Conventional amino linker and ssH amino linker



application contd ..... see back cover

# ssR Amino Linkers



## Applications of ssR linkers; ..... Contd.....

- ssH and ssMe amino linkers are stable under alkaline condition. ssH-linker is stable even in carbonate buffer.
- For solid phase labeling, transacetylation reactions were suppressed by a) keeping the terminal amino groups protonated and b) by activating the exogenous molecules before the coupling reaction.
- pKa of the aminoethyl carbamate structure is lower than the aliphatic amine, this might be responsible for the **rapid MMTr removal and efficient labeling reaction**.
- Attachment of oligonucleotides to **chips (microarray: high quality amino-ODNs)**. Immobilization of the modified oligonucleotides to various surfaces for oligonucleotide labeling with corresponding reporter molecules.

## Reaction conditions of the PEGYLATION:

### Conjugation condition with Branched 40 K PEG\*:

Amino-modified oligos of 25-bases (5  $\mu$ mole) were incubated with activated PEG-NHS\* (500  $\mu$ mole) in 250 mM phosphate buffer (pH8) at 40 °C. After 20 min, aliquots (30  $\mu$ L) removed from the solvent, was combined with distilled water, followed by HPLC analysis using anion-exchange DEAE column.

\*PEG-NHS: Shown in the example in figure is symmetrical branched PEG-NHS Ester-MW 40,000

## ssMe-linker:

1. ssMe is superior to ssH in both conjugation efficiencies and MMTr-removal.
2. The cost of ssMe-linker is higher due to single isomer (S) being used in the synthesis.
3. Degradation of ssMe-linker under alkaline conditions is slightly more as compared to the ssH linker oligos.

**Important Note:** It is preferable not to deprotect MMTr during oligo synthesis, as some carbamate cleavage during aq. ammonia deprotection can occur.

## Licensing & Trademarks

ChemGenes Corp. holds worldwide license, US7491857B2; Y. Komatsu; N. Kojima; K. Sata; Ken Nonaka; Y. Fujiwara. NIAIST & DNA chip Research Inc, Japan.

## References:

1. Komatsu, Y. et. al. *Bioorg Med. Chem.* **2008**, 16, 941.
2. Kojima, N. et. al. *Bioorg. Med. Chem. Lett.* **2009**, 19, 2144.
3. Komatsu, Y. et. al. US 2008/0227968 A1.